


Report Title

Test Name:	HLAB27-Q	Test Sample:	HLApc
Test Time:	2025-09-20-10-20-23_1	Report Time:	2025-09-20

Well	Dyes	Targets	Type	Ct	Copies	Result-1	Result-2	Final Result
Well A	FAM	HLA	Unknown	15.38	NA	HLA Positive	HLA Positive, IC Detected	HLA Positive, IC Detected
	CY5	IC	Unknown	13.38	NA	IC Detected		



When the Ct value of the sample to be tested is outside the reference range and there is a typical S-shaped amplification curve, the test result is positive;
When the Ct value of the sample to be tested is within the reference range, or the Ct value is outside the reference range but there is no typical S-shaped amplification curve, the test result is negative.