



# Other Infectious Diseases

For various viral infections, high sensitivity and specificity are required for early detection and prevention of further spread. Qualitative detection of infectious diseases is accomplished by viral RNA detection using reverse transcription–polymerase chain reaction (RT-PCR) kits.



## Sample Source

Blood/ Serum/ Plasma, Urine, Prostate fluid and Pap smear



## Target Specific

Individual marker specific kits with specific designed primers and fluorescent probes



## Easy & Flexible

Highly compatible with recent available RT-PCR instruments

### EBV-Q Real Time PCR Kit

Quantitative detection of Epstein Barr Virus 1 & 2 specific EBNA1 gene in human blood samples

### Nipah-Q Real Time PCR Kit

Quantitative detection of RNA from Nipah virus (NiV) in human specimens

### Human Parvovirus (B19)-Q Real-Time PCR Kit

Quantitative detection of Human Parvovirus (B19) specific NS1 gene in clinical sample

### Rota-Q Real-Time PCR Kit

Detection of Rotavirus A specific RNA in clinical samples

## H. Pylori-Q Real Time PCR Kit

Quantitative detection of Helicobacter pylori bacteria specific rpoD gene in clinical samples

## BK-Q Real Time PCR Kit

Quantitative detection of Polyomavirus BK virus specific VP1 gene in clinical samples

## VZV-Q Real Time PCR Kit

Quantitative detection of Varicella zoster virus (VZV) specific ORF62 gene in clinical samples

## C. Albicans (CA)-Q Real Time PCR Kit

Quantitative detection of Candida albicans specific ribonuclease P RNA (RPR1) gene in clinical samples

## EEC-Q Real Time PCR Kit

Detection of Enterovirus, Enterovirus 71 and Coxsackie Virus A16 in clinical samples including nasal and pharyngeal secretions, etc

## SS-Q Real Time PCR Kit

Detection of Salmonella & Shigella in stool/ water samples by Real Time PCR systems

Commercial Name	Cat No.	Compatible Instruments
H. Pylori-Q Real-Time PCR Kit	G2M706321	Instruments with FAM, ROX, Cy5, HEX/VIC Channels including ABI Prism®7500 & Quantstudio, BioRad CFX96, Roche Lightcycler, Qiagen Rotor-Gene, etc.
BK-Q Real Time PCR Kit	G2M705421	
EBV-Q Real Time PCR Kit	G2M705621	
Nipah-Q Real Time PCR Kit	G2M705921	
VZV-Q Real Time PCR Kit	G2M705821	
Human Parvovirus-Q Real Time PCR Kit	G2M706521	
C. Albicans-Q Real Time PCR Kit	G2M708021	
Rota-Q Real Time PCR Kit	G2M801521	
SS-Q Real Time PCR Kit	G2M801621	
EEC-Q Real Time PCR Kit	G2M801221	



### Genes2Me Pvt. Ltd.

1105, 11th Floor, Tower B, SAS Tower, Medicity,  
Sector - 38, Gurgaon - 122 001, Haryana, India  
Tel : + 91 18001 214030 / + 91 88000 23600  
E-mail : contact@genes2me.com www.genes2me.com

