

Reproductive Health

Genes2Me RT-PCR Kits offer efficient qualitative clinical diagnosis of patient specimens associated with inflammatory diseases, bacterial infection in reproductive organs, infertility and sexually transmitted infections (STIs).



Multiple Specimen

Wide variety sample specimens including Urine, Tissue Swab (Vagina/ cervix/ urethra), Endocervical Swab, Genital Swab/ Fluid and Prostatic Fluid



High Specificity

Corresponding better to quality assurance including Internal Control (IC) and no cross reactivity with other pathogens

Syphilis-Q Real Time PCR Kit

Qualitative detection of Treponema pallidum specific gene with Taqman Chemistry on Real Time PCR based methodology

HSV-Q Real Time PCR Kit

Quantitative detection of Herpes Simplex Virus 1 (HSV-1) and Herpes Simplex Virus 2 (HSV-2) specific UL30 gene in clinical samples

STI7-Q Real Time PCR Kit

Qualitative detection of seven major sexually transmitted infections (STIs) causative pathogens in a single Real-Time PCR reaction from clinical samples

N. Gonorrhoeae (NG)-Q Real Time PCR Kit

Quantitative detection of Neisseria Gonorrhoeae strain specific Type I restriction modification DNA specificity domain-containing protein in sample sources

T. Vaginalis (TV)-Q Real Time PCR Kit

Quantitative detection of Trichomonas vaginalis (TV) specific PMS1 gene in human samples

M. Hominis (MH)-Q Real Time PCR Kit

Quantitative detection of Mycoplasma hominis (MH) bacterial strain specific gap gene in clinical samples

U. Urealyticum-Q Real Time PCR Kit

In vitro quantification of Ureaplasma urealyticum in clinical samples using target sequence within the urease complex

HHV-Q Real Time PCR Kit

Identification and differentiation of Human Herpesvirus variants in clinical samples

Commercial Name	Cat No.	Compatible Instruments
Neisseria Gonorrhoeae (NG) Q Real Time PCR Kit	G2M706421	Instruments with FAM, ROX, Cy5, HEX/VIC Channels including ABI Prism@7500 & Quantstudio, BioRad CFX96, Roche Lightcycler, Qiagen Rotor-Gene, etc.
Trichomonas Vaginalis (TV) Q Real Time PCR Kit	G2M706621	
Mycoplasma Hominis (MH) Q Real Time PCR Kit	G2M707221	
Syphilis-Q Real Time PCR Kit	G2M705521	
HSV-Q Real Time PCR Test Kit	G2M705221	
STI7-Q Real Time PCR Kit	G2M709121	
U. Urealyticum-Q Real Time PCR Kit	G2M801321	
HHV-Q Real Time PCR Kit	G2M802321	



Genes2Me Pvt. Ltd.

1105, 11th Floor, Tower B, SAS Tower, Medicity,
Sector - 38, Gurgaon - 122 001, Haryana, India
Tel : + 91 18001 214030 / + 91 88000 23600
E-mail : contact@genes2me.com www.genes2me.com

